



KLH31 mouse mAb

Catalog No	YP-mAb-10066
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KLHL31 BKLHD6 KBTBD1 KLHL
Protein Name	Kelch-like protein 31 (BTB and kelch domain-containing protein 6) (Kelch repeat and BTB domain-containing protein 1) (Kelch-like protein KLHL)
Immunogen	Synthesized peptide derived from human KLH31 AA range: 342-392
Specificity	This antibody detects endogenous levels of human KLH31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Strongly expressed in skeletal muscle and weakly in heart. According to PubMed:15302408, not expressed in other tissues. According to PubMed:18719355, abundantly expressed in both embryonic skeletal and heart tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images