



ME3L2 mouse mAb

Catalog No	YP-mAb-10109
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MEIS3P2
Protein Name	Putative homeobox protein Meis3-like 2
Immunogen	Synthesized peptide derived from human ME3L2 AA range: 200-250
Specificity	This antibody detects endogenous levels of human ME3L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images