



FBX31 mouse mAb

Catalog No	YP-mAb-10115
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FBXO31 FBX14 FBX31 PP2386
Protein Name	F-box only protein 31
Immunogen	Synthesized peptide derived from human FBX31 AA range: 17-67
Specificity	This antibody detects endogenous levels of human FBX31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Highly expressed in brain. Expressed at moderate levels in most tissues, except bone marrow.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images