



PDE6C mouse mAb

Catalog No	YP-mAb-10122
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PDE6C PDEA2
Protein Name	Cone cGMP-specific 3',5'-cyclic phosphodiesterase subunit alpha' (EC 3.1.4.35) (cGMP phosphodiesterase 6C)
Immunogen	Synthesized peptide derived from human PDE6C AA range: 474-524
Specificity	This antibody detects endogenous levels of human PDE6C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images