



ADRO mouse mAb

Catalog No	YP-mAb-10180
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	FDXR ADXR
Protein Name	NADPH:adrenodoxin oxidoreductase, mitochondrial (AR) (Adrenodoxin reductase) (EC 1.18.1.2) (Ferredoxin--NADP(+) reductase) (Ferredoxin reductase)
Immunogen	Synthesized peptide derived from human ADRO AA range: 364-414
Specificity	This antibody detects endogenous levels of human ADRO
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	54kD
Cell Pathway	Mitochondrion inner membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images