



TT30B mouse mAb

Catalog No	YP-mAb-10192
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TTC30B
Protein Name	Tetratricopeptide repeat protein 30B (TPR repeat protein 30B)
Immunogen	Synthesized peptide derived from human TT30B AA range: 423-473
Specificity	This antibody detects endogenous levels of human TT30B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell projection, cilium .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images