



AP3B1 mouse mAb

Catalog No	YP-mAb-10194
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AP3B1 ADTB3A
Protein Name	AP-3 complex subunit beta-1 (Adapter-related protein complex 3 subunit beta-1) (Adaptor protein complex AP-3 subunit beta-1) (Beta-3A-adaptin) (Clathrin assembly protein complex 3 beta-1 large chain)
Immunogen	Synthesized peptide derived from human AP3B1 AA range: 525-575
Specificity	This antibody detects endogenous levels of human AP3B1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, clathrin-coated vesicle membrane ; Peripheral membrane protein ; Cytoplasmic side . Golgi apparatus . Component of the coat surrounding the cytoplasmic face of coated vesicles located at the Golgi complex. .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images