



WIT1 Polyclonal Antibody

Catalog No	YP-mAb-10898
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	WT1-AS WIT1
Protein Name	WIT1
Immunogen	Synthesized peptide derived from human WIT1 AA range: 20-100
Specificity	This antibody detects endogenous levels of human WIT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Putative Wilms tumor upstream neighbor 1 gene protein (WIT-1;Wilms tumor-associated antisense RNA)
Observed Band	
Cell Pathway	
Tissue Specificity	Kidney,Placenta,Uterus,
Function	cell proliferation,
Background	This gene is located upstream of the Wilms tumor 1 (WT1) gene; these two genes are bi-directionally transcribed from the same promoter region. This gene is imprinted in kidney, with preferential expression from the paternal allele. Imprinting defects at chromosome 11p13 may contribute to tumorigenesis. [provided by RefSeq, May 2014],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images