



IGF1 mouse mAb

Catalog No	YP-mAb-12454
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	IGF1 IBP1
Protein Name	IGF1
Immunogen	Synthesized peptide derived from human IGF1
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat IGF1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Insulin-like growth factor I (IGF-I;Mechano growth factor;MGF;Somatomedin-C)
Observed Band	21kD
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images