



IL33 mouse mAb

Catalog No	YP-mAb-12460
Isotype	IgG
Reactivity	Mouse
Applications	WB
Gene Name	IL33
Protein Name	IL33
Immunogen	Synthesized peptide derived from mouse IL33
Specificity	This antibody detects endogenous levels of Mouse IL33
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Interleukin-33 (IL-33;Interleukin-1 family member 11;IL-1F11;Nuclear factor from high endothelial venules;NF-HEV)
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images