



IL27B mouse mAb

Catalog No	YP-mAb-12468
Isotype	IgG
Reactivity	Mouse
Applications	WB
Gene Name	IL27B
Protein Name	IL27B
Immunogen	Synthesized peptide derived from mouse IL27B
Specificity	This antibody detects endogenous levels of Mouse IL27B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Interleukin-27 subunit beta (IL-27 subunit beta;IL-27B;Epstein-Barr virus-induced gene 3 protein;EBV-induced gene 3 protein)
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images