



SNAT mouse mAb

Catalog No	YP-mAb-12511
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	AANAT SNAT
Protein Name	SNAT
Immunogen	Synthesized peptide derived from human SNAT
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat SNAT
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Serotonin N-acetyltransferase (Serotonin acetylase;EC 2.3.1.87;Aralkylamine N-acetyltransferase;AA-NAT)
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Highly expressed in pineal gland and at lower levels in the retina. Weak expression in several brain regions and in the pituitary gland.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images