



IRAK4 mouse mAb

Catalog No	YP-mAb-12532
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	IRAK4
Protein Name	IRAK4
Immunogen	Synthesized peptide derived from human IRAK4
Specificity	This antibody detects endogenous levels of IRAK4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.14% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Interleukin-1 receptor-associated kinase 4 (IRAK-4) (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-64)
Observed Band	55KD
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images