



HIVEP2 mouse mAb

Catalog No	YP-mAb-12563
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	HIVEP2
Protein Name	HIVEP2
Immunogen	Synthesized peptide derived from human HIVEP2
Specificity	This antibody detects endogenous levels of HIVEP2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.54% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Transcription factor HIVEP2 (Human immunodeficiency virus type I enhancer-binding protein 2) (HIV-EP2) (MHC-binding protein 2) (MBP-2)
Observed Band	
Cell Pathway	Nucleus.
Tissue Specificity	Expressed in brain and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images