



# sesn2 mouse mAb

|                                  |   |
|----------------------------------|---|
| <b>Catalog No</b>                | YP-mAb-12582  |
| <b>Isotype</b>                   | IgG   |
| <b>Reactivity</b>                | Human; Mouse;Rat  |
| <b>Applications</b>              | WB  |
| <b>Gene Name</b>                 | SESN2 SEST2   |
| <b>Protein Name</b>              | sesn2   |
| <b>Immunogen</b>                 | Synthesized peptide derived from human sesn2  |
| <b>Specificity</b>               | This antibody detects endogenous levels of sesn2 at Human, Mouse,Rat  |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, and 0.74% sodium azide.  |
| <b>Source</b>                    | Monoclonal, Mouse,IgG   |
| <b>Purification</b>              | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.              |
| <b>Dilution</b>                  | WB 1:500-1:2000   |
| <b>Concentration</b>             | 1 mg/ml   |
| <b>Purity</b>                    | ≥90%  |
| <b>Storage Stability</b>         | -20°C/1 year  |
| <b>Synonyms</b>                  | Sestrin-2 (Hi95)  |
| <b>Observed Band</b>             |   |
| <b>Cell Pathway</b>              | Cytoplasm .   |
| <b>Tissue Specificity</b>        | Widely expressed.   |
| <b>Function</b>                  |   |
| <b>Background</b>                |   |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!  |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images