



G-CSFR Monoclonal Antibody

Catalog No	YP-mAb-13714
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CSF3R
Protein Name	Granulocyte colony-stimulating factor receptor
Immunogen	Synthesized peptide derived from the Internal region of human G-CSFR.
Specificity	G-CSFR Monoclonal Antibody detects endogenous levels of G-CSFR protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CSF3R; GCSFR; Granulocyte colony-stimulating factor receptor; G-CSF receptor; G-CSF-R; CD antigen CD114
Observed Band	92156DA
Cell Pathway	[Isoform 2]: Secreted .; Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	One or several isoforms have been found in myelogenous leukemia cell line KG-1, leukemia U-937 cell line, in bone marrow cells, placenta, and peripheral blood granulocytes. Isoform GCSFR-2 is found only in leukemia U-937 cells. Isoform GCSFR-3 is highly expressed in placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images