



MAPKAPK-2 (phospho-Thr334) mouse mAb

Catalog No	YP-mAb-14599
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MAPKAPK2
Protein Name	MAPKAPK-2 (Thr334)
Immunogen	Synthesized phospho peptide around human MAPKAPK-2 (Thr334)
Specificity	This antibody detects endogenous levels of Human Mouse Rat MAPKAPK-2 (phospho-Thr334)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAP kinase-activated protein kinase 2 (MAPK-activated protein kinase 2) (MAPKAP kinase 2) (MAPKAP-K2) (MAPKAPK-2) (MK-2) (MK2) (EC 2.7.11.1)
Observed Band	57305DA
Cell Pathway	Cytoplasm . Nucleus . Phosphorylation and subsequent activation releases the autoinhibitory helix, resulting in the export from the nucleus into the cytoplasm.
Tissue Specificity	Expressed in all tissues examined.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images