



PCYT2 mouse mAb

Catalog No	YP-mAb-17123
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	type（类型）： Monoclonal Antibody
Protein Name	
Immunogen	Recombinant fusion protein of human PCYT2 (NP_002852.1).
Specificity	matters needing attention（注意事项）： Avoid repeated freezing and thawing!
Formulation	
Source	Host（宿主）： Mouse
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	Detects a band of approximately 44kD
Cell Pathway	
Tissue Specificity	Usage suggestions（使用建议）： This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.
Function	Disclaimer note（免责声明）： The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images