



VEGF Receptor 1 mouse mAb

Catalog No	YP-mAb-17140
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	FLT1
Protein Name	
Immunogen	A synthesized peptide derived from human VEGF Receptor 1
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Modification（修改）Unmodified
Observed Band	Form（类型）Liquid
Cell Pathway	
Tissue Specificity	Buffer System（缓冲系统） Mouse IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Function	
Background	Gene ID: 2321
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images