



Collagen II α 1 (Cleaved-Ala1241) mouse mAb

Catalog No	YP-mAb-17141
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	COL2A1
Protein Name	Collagen II α 1 (Cleaved-Ala1241) Human Gene Id: 1280 Human Swiss Prot No: P02458 Mouse Gene Id: 12824 Mouse Swiss Prot No: P28481 Rat Gene Id: 25412 Rat Swiss Prot No: P05539
Immunogen	Synthesized peptide derived from human Collagen II α 1 (Cleaved-Ala1241)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat Collagen II α 1 (Cleaved-Ala1241, protein was cleaved amino acid sequence between 1241-1242)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	

**Cell Pathway****Tissue Specificity****Function****Background****matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images