



ERAL1 mouse mAb

Catalog No	YP-mAb-17236
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	WB
Gene Name	ERAL1;HERA
Protein Name	GTPase Era, mitochondrial;Conserved ERA-like GTPase (CEGA)
Immunogen	Synthesized peptide derived from human C-terminal ERAL1
Specificity	This antibody detects endogenous levels of ERAL1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	#N/A
Observed Band	48kD
Cell Pathway	#N/A
Tissue Specificity	#N/A
Function	#N/A
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images