



TUTLA mouse mAb

Catalog No	YP-mAb-17281
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	WB
Gene Name	IGSF9;IGSF9A; KIAA1355; NRT1
Protein Name	Protein turtle homolog A;Immunoglobulin superfamily member 9A (IgSF9A)
Immunogen	Synthesized peptide derived from human C-ternal TUTLA
Specificity	This antibody detects endogenous levels of TUTLA at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Mouse,Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	#N/A
Observed Band	130kD
Cell Pathway	#N/A
Tissue Specificity	#N/A
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images