



EPHX2 mouse mAb

Catalog No	YP-mAb-17670
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	WB
Gene Name	EPHX2
Protein Name	
Immunogen	
Specificity	EPHX2
Formulation	<p>Bifunctional epoxide hydrolase 2 [Includes: Cytosolic epoxide hydrolase 2 (CEH) (EC 3.3.2.10) (Epoxide hydratase) (Soluble epoxide hydrolase) (SEH); Lipid-phosphate phosphatase (EC 3.1.3.76)]</p> <p>Human Gene Id: 2053</p> <p>Human Swiss Prot No: P34913</p> <p>Mouse Gene Id: 13850</p> <p>Mouse Swiss Prot No: P34914</p> <p>Rat Gene Id: 65030</p> <p>Rat Swiss Prot No: P80299</p>
Source	Mouse, Monoclonal
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	The antibody was affinity-purified from Mouse serum by affinity-chromatography using specific immunogen.

**Observed Band****Cell Pathway****Tissue Specificity** 61kD**Function****Background****matters needing attention** Avoid repeated freezing and thawing!**Usage suggestions** This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.**Products Images**