



BDH2 mouse mAb

Catalog No	YP-mAb-17720
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	BDH2
Protein Name	
Immunogen	Full length fusion protein
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DHRS6; EFA6R; SDR15C1; UCPA-OR; UNQ6308; PRO20933
Observed Band	Calculated MW: 27 kDa; Observed MW: 27 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Dehydrogenase that mediates the formation of 2,5-dihydroxybenzoic acid (2,5-DHBA), a siderophore that shares structural similarities with bacterial enterobactin and associates with LCN2, thereby playing a key role in iron assimilation and homeostasis.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images