



Phospho-Tau (Ser356) mouse mAb

Catalog No	YP-mAb-17775
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	MAPT
Protein Name	
Immunogen	Synthetic peptide of human MAPT
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	Calculated MW: 79 kDa; Observed MW: 48,62,78 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	This gene encodes the microtubule-associated protein tau (MAPT) whose transcript undergoes complex, regulated alternative splicing, giving rise to several mRNA species. MAPT transcripts are differentially expressed in the nervous system, depending on stage of neuronal maturation and neuron type.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images