



SET mouse mAb

Catalog No	YP-mAb-17842
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SET
Protein Name	
Immunogen	A synthesized peptide derived from human SET
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	2PP2A; IGAAD; TAF-I; I2PP2A; IPP2A2; PHAPII; TAF-IBETA
Observed Band	Calculated MW: 33 kDa; Observed MW: 39 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Multitasking protein, involved in apoptosis, transcription, nucleosome assembly and histone binding. Isoform 2 anti-apoptotic activity is mediated by inhibition of the GZMA-activated DNase, NME1.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images