



CCR6 mouse mAb

Catalog No	YP-mAb-17863
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	CCR6
Protein Name	
Immunogen	Synthetic peptide of human CCR6
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCR6; CKRL3; CMKBR6; GPR29; STRL22; C-C chemokine receptor type 6; C-C CKR-6; CC-CKR-6; CCR-6; Chemokine receptor-like 3; CKR-L3; DRY6; G-protein coupled receptor 29; GPR-CY4; GPRCY4; LARC receptor; CD196
Observed Band	Calculated MW: 42 kDa; Observed MW: 42 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Receptor for the C-C type chemokine CCL20.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images