



LGI2 mouse mAb

Catalog No	YP-mAb-17936
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LGI2 KIAA1916 LGIL2
Protein Name	Leucine-rich repeat LGI family member 2 (LGI1-like protein 2) (Leucine-rich glioma-inactivated protein 2)
Immunogen	Synthesized peptide derived from human LGI2
Specificity	This antibody detects endogenous levels of LGI2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	60kD
Cell Pathway	Secreted .
Tissue Specificity	Brain, heart and placenta.
Function	Required for the development of soma-targeting inhibitory GABAergic synapses made by parvalbumin-positive basket cells.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images