



MIOX mouse mAb

Catalog No	YP-mAb-17944
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MIOX ALDRL6 KSP32 RSOR
Protein Name	Inositol oxygenase (EC 1.13.99.1) (Aldehyde reductase-like 6) (Kidney-specific protein 32) (Myo-inositol oxygenase) (MI oxygenase) (Renal-specific oxidoreductase)
Immunogen	Synthesized peptide derived from human MIOX
Specificity	This antibody detects endogenous levels of MIOX at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	31kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Kidney specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images