



SSX4 mouse mAb

Catalog No	YP-mAb-18020
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	SSX4 SSX4A; SSX4B
Protein Name	Protein SSX4 (Cancer/testis antigen 5.4) (CT5.4)
Immunogen	Synthesized peptide derived from human SSX4
Specificity	This antibody detects endogenous levels of SSX4 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	21kD
Cell Pathway	
Tissue Specificity	
Function	Could act as a modulator of transcription.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images