



SPANXB1 mouse mAb

Catalog No	YP-mAb-18021
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	SPANXB1 SPANXB; SPANXB2; SPANXF1
Protein Name	Sperm protein associated with the nucleus on the X chromosome B/F (Cancer/testis antigen 11.2) (CT11.2) (Nuclear-associated protein SPAN-Xb) (SPANX-B) (Nuclear-associated protein SPAN-Xf) (SPANX-F)
Immunogen	Synthesized peptide derived from human SPANXB1
Specificity	This antibody detects endogenous levels of SPANXB1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	11kD
Cell Pathway	Cytoplasm . Nucleus . Associated with nuclear craters.
Tissue Specificity	Detected in testis and sperm.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images