



SCUBE3 mouse mAb

Catalog No	YP-mAb-18086
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SCUBE3 CEGF3
Protein Name	Signal peptide, CUB and EGF-like domain-containing protein 3
Immunogen	Synthesized peptide derived from human SCUBE3
Specificity	This antibody detects endogenous levels of SCUBE3 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	109kD
Cell Pathway	Secreted . Cell surface .
Tissue Specificity	Highly expressed in osteoblasts. In normal lung, mainly expressed in bronchial epithelial cells. Tends to be up-regulated in lung cancer cells.
Function	Is a positive regulator of the BMP signaling pathway, required for proper chondrogenesis, osteogenesis and skeletal development. It acts as coreceptor for BMP ligands, particularly BMP2 and BMP4, facilitating their interactions with BMP type I receptors . It is required for ligand-induced recruitment of BMP receptors to lipid rafts (By similarity). Binds to TGFB2 and activates TGFB signaling. In lung cancer cells, could serve as an endogenous autocrine and paracrine ligand of TGFB2, which could regulate TGFB2 signaling and hence modulate epithelial-mesenchymal transition and cancer progression.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images