



TRIM32 mouse mAb

Catalog No YP-mAb-18101

Isotype IgG

Reactivity Human;Mouse

Applications WB

Gene Name

Protein Name

Immunogen

Specificity

Formulation

Source

Purification The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.

Dilution WB 1:500-1:2000

Concentration 1 mg/ml

Purity ≥90%

Storage Stability -20°C/1 year

Synonyms

Observed Band

Cell Pathway

Tissue Specificity

Function

Background

matters needing attention Avoid repeated freezing and thawing!

Usage suggestions This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images