



COMMD8 mouse mAb

Catalog No	YP-mAb-18171
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	COMMD8 MDS022
Protein Name	COMM domain-containing protein 8
Immunogen	Synthesized peptide derived from human COMMD8
Specificity	This antibody detects endogenous levels of COMMD8 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	20kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Widely expressed with highest expression in thyroid.
Function	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes . May down-regulate activation of NF-kappa-B .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images