



BHLHA15 mouse mAb

Catalog No	YP-mAb-18196
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	BHLHA15 BHLHB8 MIST1
Protein Name	Class A basic helix-loop-helix protein 15 (bHLHa15) (Class B basic helix-loop-helix protein 8) (bHLHb8) (Muscle, intestine and stomach expression 1) (MIST-1)
Immunogen	Synthesized peptide derived from human BHLHA15
Specificity	This antibody detects endogenous levels of BHLHA15 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	21kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in brain, liver, spleen and skeletal muscle.
Function	Plays a role in controlling the transcriptional activity of MYOD1, ensuring that expanding myoblast populations remain undifferentiated. Repression may occur through muscle-specific E-box occupancy by homodimers. May also negatively regulate bHLH-mediated transcription through an N-terminal repressor domain. Serves as a key regulator of acinar cell function, stability, and identity. Also required for normal organelle localization in exocrine cells and for mitochondrial calcium ion transport. May function as a unique regulator of gene expression in several different embryonic and postnatal cell lineages. Binds to the E-box consensus sequence 5'-CANNTG-3' (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images