




LRP5L mouse mAb

Catalog No	YP-mAb-11347
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LRP5L
Protein Name	LRP5L
Immunogen	Synthesized peptide derived from human LRP5L AA range: 68-118
Specificity	This antibody detects endogenous levels of LRP5L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	plasma membrane,receptor complex,
Tissue Specificity	
Function	similarity:Contains 5 LDL-receptor class B repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using LRP5L mouse mAb