



# HYAL1 Rabbit pAb

<b>Catalog No</b>	YP-Ab-18569
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Mouse,Rat
<b>Applications</b>	WB
<b>Gene Name</b>	HYAL1 LUCA1
<b>Protein Name</b>	Hyaluronidase-1 (Hyal-1) (Hyaluronoglucosaminidase-1) (Lung carcinoma protein 1) (LuCa-1)
<b>Immunogen</b>	Synthesized peptide derived from human HYAL1
<b>Specificity</b>	This antibody detects endogenous levels of HYAL1 at Human, Mouse,Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	48kD
<b>Cell Pathway</b>	
<b>Tissue Specificity</b>	
<b>Function</b>	May have a role in promoting tumor progression. May block the TGFB1-enhanced cell growth.
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

