



HYAL2 Rabbit pAb

Catalog No	YP-Ab-18570
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HYAL1 LUCA1
Protein Name	Hyaluronidase-1 (Hyal-1) (Hyaluronoglucosaminidase-1) (Lung carcinoma protein 1) (LuCa-1)
Immunogen	Synthesized peptide derived from human HYAL2
Specificity	This antibody detects endogenous levels of HYAL2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	52kD
Cell Pathway	
Tissue Specificity	
Function	Hydrolyzes high molecular weight hyaluronic acid to produce an intermediate-sized product which is further hydrolyzed by sperm hyaluronidase to give small oligosaccharides. Displays very low levels of activity. Associates with and negatively regulates MST1R.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images