




Cdc50C Polyclonal Antibody

Catalog No	YP-Ab-16688
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IF;ELISA
Gene Name	TMEM30C
Protein Name	Cell cycle control protein 50C
Immunogen	The antiserum was produced against synthesized peptide derived from human TMEM30C. AA range:51-100
Specificity	Cdc50C Polyclonal Antibody detects endogenous levels of Cdc50C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMEM30C; CDC50C; Cell cycle control protein 50C; Transmembrane protein 30C
Observed Band	
Cell Pathway	Membrane ; Single-pass membrane protein .
Tissue Specificity	Specifically expressed in testis.
Function	
Background	TMEM30C, also named as CDC50C, belongs to the CDC50/LEM3 family. TMEM30C is specifically expressed in testis.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunofluorescence analysis of MCF7 cells, using TMEM30C Antibody. The picture on the right is blocked with the synthesized peptide.