



PPM1G rabbit pAb

Catalog No	YP-Ab-17234
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	PPM1G PPM1C
Protein Name	Protein phosphatase 1G (EC 3.1.3.16) (Protein phosphatase 1C) (Protein phosphatase 2C isoform gamma) (PP2C-gamma) (Protein phosphatase magnesium-dependent 1 gamma)
Immunogen	Synthesized peptide derived from human C-terminal PPM1G
Specificity	This antibody detects endogenous levels of PPM1G at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein phosphatase 1G (EC 3.1.3.16) (Protein phosphatase 1C) (Protein phosphatase 2C isoform gamma) (PP2C-gamma) (Protein phosphatase magnesium-dependent 1 gamma)
Observed Band	
Cell Pathway	Cytoplasm . Membrane ; Lipid-anchor .
Tissue Specificity	Widely expressed. Most abundant in testis, skeletal muscle, and heart.
Function	0
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images