



ERAL1 rabbit pAb

Catalog No	YP-Ab-17236
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	IHC,WB
Gene Name	ERAL1;HERA
Protein Name	GTPase Era, mitochondrial;Conserved ERA-like GTPase (CEGA)
Immunogen	Synthesized peptide derived from human C-ternal ERAL1
Specificity	This antibody detects endogenous levels of ERAL1 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit,polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ERROR 42
Observed Band	48kD
Cell Pathway	ERROR 42
Tissue Specificity	ERROR 42
Function	ERROR 42
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

