



G6PC rabbit pAb

Catalog No	YP-Ab-17245
Isotype	IgG
Reactivity	Human, Mouse,Rat
Applications	IHC,WB
Gene Name	G6PC G6PT
Protein Name	Glucose-6-phosphatase (G-6-Pase) (G6Pase) (EC 3.1.3.9) (Glucose-6-phosphatase alpha) (G6Pase-alpha)
Immunogen	Synthesized peptide derived from human C-ternal G6PC
Specificity	This antibody detects endogenous levels of G6PC at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit,polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glucose-6-phosphatase (G-6-Pase) (G6Pase) (EC 3.1.3.9) (Glucose-6-phosphatase alpha) (G6Pase-alpha)
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	Hydrolyzes glucose-6-phosphate to glucose in the endoplasmic reticulum. Forms with the glucose-6-phosphate transporter (SLC37A4/G6PT) the complex responsible for glucose production in the terminal step of glycogenolysis and gluconeogenesis. Hence, it is the key enzyme in homeostatic regulation of blood glucose levels.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images