



TFAP4 rabbit pAb

Catalog No	YP-Ab-17249
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	TFAP4 BHLHC41
Protein Name	Transcription factor AP-4 (Activating enhancer-binding protein 4) (Class C basic helix-loop-helix protein 41) (bHLHc41)
Immunogen	Synthesized peptide derived from human N-terminal TFAP4
Specificity	This antibody detects endogenous levels of TFAP4 at Human
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Transcription factor AP-4 (Activating enhancer-binding protein 4) (Class C basic helix-loop-helix protein 41) (bHLHc41)
Observed Band	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	Transcription factor that activates both viral and cellular genes by binding to the symmetrical DNA sequence 5'-CAGCTG-3'.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

