



ITIH5 rabbit pAb

Catalog No	YP-Ab-17260
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	ITIH5 KIAA1953 PP14776 UNQ311/PRO354
Protein Name	Inter-alpha-trypsin inhibitor heavy chain H5 (ITI heavy chain H5) (ITI-HC5) (Inter-alpha-inhibitor heavy chain 5)
Immunogen	Synthesized peptide derived from human C-terminal ITIH5
Specificity	This antibody detects endogenous levels of ITIH5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Inter-alpha-trypsin inhibitor heavy chain H5 (ITI heavy chain H5) (ITI-HC5) (Inter-alpha-inhibitor heavy chain 5)
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	Abundantly expressed in placenta. Less abundant expression in mammary gland and ovary. Expression is barely detectable levels in all other tissues tested.
Function	May act as a tumor suppressor.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

