



METL8 rabbit pAb

Catalog No	YP-Ab-17278
Isotype	IgG
Reactivity	Human, Mouse
Applications	IHC, WB
Gene Name	METTL8
Protein Name	Methyltransferase-like protein 8 (EC 2.1.1.-)
Immunogen	Synthesized peptide derived from human C-terminal METL8
Specificity	This antibody detects endogenous levels of METL8 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Rabbit, polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000 IHC 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Methyltransferase-like protein 8 (EC 2.1.1.-)
Observed Band	
Cell Pathway	Mitochondrion . Mitochondrial protein: the cytoplasmic or nuclear localization observed by some groups is either the result of an incorrect localization caused by N-terminal tagging that interferes with mitochondrial targeting, or splice isoforms that lack the N-terminal mitochondrial transit sequence. .
Tissue Specificity	
Function	Mitochondrial S-adenosyl-L-methionine-dependent methyltransferase that mediates N(3)-methylcytidine modification of residue 32 of the tRNA anticodon loop of mitochondrial tRNA(Ser)(UCN) and tRNA(Thr) . N(3)-methylcytidine methylation modification regulates mitochondrial translation efficiency and is required for activity of the respiratory chain . N(3)-methylcytidine methylation of mitochondrial tRNA(Ser)(UCN) requires the formation of N(6)-dimethylallyl-adenosine(37) (i6A37) by TRIT1 as prerequisite . May also mediate N(3)-methylcytidine modification of mRNAs . The existence of N(3)-methylcytidine modification on mRNAs is however unclear, and additional evidences are required to confirm the role of the N(3)-methylcytidine-specific mRNA methyltransferase activity of METTL8 in vivo .
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images