



# TUTLA rabbit pAb

<b>Catalog No</b>	YP-Ab-17281
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human, Mouse,Rat
<b>Applications</b>	IHC,WB
<b>Gene Name</b>	IGSF9;IGSF9A; KIAA1355; NRT1
<b>Protein Name</b>	Protein turtle homolog A;Immunoglobulin superfamily member 9A (IgSF9A)
<b>Immunogen</b>	Synthesized peptide derived from human C-ternal TUTLA
<b>Specificity</b>	This antibody detects endogenous levels of TUTLA at Human, Mouse,Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Rabbit,polyclonal
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Dilution</b>	WB 1:500-2000 IHC 1:50-200
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	ERROR 42
<b>Observed Band</b>	130kD
<b>Cell Pathway</b>	ERROR 42
<b>Tissue Specificity</b>	ERROR 42
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

