



RHBDF2 rabbit pAb

Catalog No	YP-Ab-17691
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RHBDF2 IRHOM2 RHBDL5 RHBDL6
Protein Name	Inactive rhomboid protein 2 (iRhom2) (Rhomboid 5 homolog 2) (Rhomboid family member 2) (Rhomboid veinlet-like protein 5) (Rhomboid veinlet-like protein 6)
Immunogen	Synthesized peptide derived from human RHBDF2
Specificity	TThis antibody detects endogenous levels of RHBDF2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, and 0.1157% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Observed Band	94kD
Cellular localization	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Cell membrane .
Function	Regulates ADAM17 protease, a sheddase of the epidermal growth factor (EGF) receptor ligands and TNF, thereby plays a role in sleep, cell survival, proliferation, migration and inflammation. Does not exhibit any protease activity on its own.
tissue expression	Found in the epidermis and esophageal epithelium.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



