



# MFF Rabbit mAb

Catalog No	YP-Ab-17717
Isotype	IgG
Gene Name	MFF
Alternative Names	C2orf33; GL004; Mff
Research Field	Tags & Cell Markers
Product Categories	Primary antibody
Host	Rabbit
Reactivity	Human,Mouse,Rat
Application	WB
Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 34 kDa
Clonality	Monoclonal Antibody
Clonality No.	R07-5Y6
Immunogen	A synthesized peptide derived from human MFF
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Plays a role in mitochondrial and peroxisomal fission.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

