



# ARRB2 Rabbit pAb

Catalog No	YP-Ab-17718
Isotype	IgG
Gene Name	ARRB2
Alternative Names	arrestin beta 2; ARB2; ARR2; BARR2
Research Field	Signal Transduction
Product Categories	Primary antibody
Host	Rabbit
Reactivity	Human,Mouse,Rat
Application	WB
Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 46 kDa; Observed MW: 46 kDa
Clonality	Polyclonal Antibody
Immunogen	Synthetic peptide of human ARRB2
Purification	Affinity Purified
Conjugation	Unconjugated
Modification	Unmodified
Form	Liquid
Buffer System	pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Background	Functions in regulating agonist-mediated G-protein coupled receptor (GPCR) signaling by mediating both receptor desensitization and resensitization processes.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images

