



AXUD1 rabbit pAb

货号	YP-Ab-17728
同位型	IgG
靶点	CSRNP1
基因名称	CSRNP1 AXUD1 TAIP3
蛋白名称	Cysteine/serine-rich nuclear protein 1 (CSRNP-1) (Axin-1 up-regulated gene 1 protein) (Protein URAX1) (TGF-beta-induced apoptosis protein 3) (TAIP-3)
免疫原	Synthesized peptide derived from human AXUD1
特异性	This antibody detects endogenous levels of AXUD1 at Human, Mouse
组成	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
来源	Polyclonal, Rabbit,IgG
稀释	WB 1:500-2000
纯化工艺	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
分子量	65kD
功能	Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity (By similarity). May have a tumor-suppressor function. May play a role in apoptosis.
细胞定位	Nucleus
组织表达	Ubiquitous. Most abundantly expressed in lung, placenta, skeletal muscle, pancreas and leukocyte. Frequently down-regulated in lung, kidney, liver and colon cancers compared with their corresponding normal tissues.
浓度	1 mg/ml
储存	-15°C to -25°C/1 year(Do not lower than -25°C)
有关注意事项	Avoid repeated freezing and thawing!
使用建议	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images